

# Striate+™

## Instructions for use

### Implantable Collagen Membrane

Striate+™ is a biocompatible, sterile, resorbable collagen barrier membrane intended for use in guided bone and guided tissue regeneration procedures.

### Indications for use

Striate+™ is indicated for use in:

- Augmentation around implants placed in immediate extraction sockets;
- Augmentation around implants placed in delayed extraction sockets;
- Filling of bone defects after root resection or removal of retained teeth;
- Guided bone regeneration in dehiscence defects; and
- Guided tissue regeneration procedures in intra-bony periodontal defects.

### Description

Striate+™ is composed of purified collagen of porcine origin and is produced using quality-controlled manufacturing processes. Raw materials are selectively sourced from within Australia from veterinary-certified animals and manufactured without crosslinking agents. Striate+™ is packed in double blister packs and sterilized by irradiation.

### Properties

Striate+™ is a barrier membrane designed to protect the bone defect space from ingrowth of gingival tissue and provide a favorable environment for osteogenesis, and to persist long enough to allow sufficient time for bone regeneration to occur.

Striate+™ collagen membranes have a bilayer structure with a rough and a smooth side. The rough side, which is placed facing the bone defect, is composed of a loose distribution of collagen bundles that provides an open scaffold which allows entry of osteogenic cells. The smooth side, which faces the gingival tissue, is composed of parallel arrangements of densely packed collagen bundles that permit passage of fluids but act as a barrier to ingrowth of epithelial cells into the defect. Animal studies have shown that Striate+™ provides an effective barrier, allowing bone regeneration to occur in the defect space, and is then completely resorbed into the surrounding tissue.

Striate+™ retains its structural integrity when wet, while conforming to the contours of the defect. Striate+™ has sufficient tensile strength to be sutured or pinned in place if required.

### Contraindications

Striate+™ should not be used if there is evidence of active infection at the treatment site.

The use of Striate+™ in patients with known sensitivity to porcine-derived materials or collagen is contraindicated.

### Directions for use

General principles of surgical practice and aseptic technique should be adhered to.

1. Prior to guided bone and tissue regeneration procedures, anti-infective therapy to eradicate any bacterial infection and counseling of the patient in good oral hygiene is highly recommended.
2. Surgically expose the bone defect and create a mucoperiosteal flap suitable for wound closure. Debride and plane the root surface carefully. Adequate debridement and implant surface disinfection should be achieved before bone augmentation around implants in peri-implantitis bone defects.
3. Fill the bone defect with bone graft or other void-filling material, taking care not to over-fill the defect. Void fillers and implants must be adequately localized or fixed in place prior to application of the Striate+™ membrane.
4. Trim Striate+™ to the required size using sterile technique. The membrane should significantly overlap the walls of the defect to assure adequate enclosure and prevent soft-tissue invasion. Striate+™ does not require pre-wetting.
5. Apply Striate+™ over the defect and apply gentle pressure until the membrane is uniformly wet and conforming and adhering to the underlying surface. The rough side of the membrane is placed facing the bone defect and the smooth side faces the oral cavity.
6. To avoid the formation of excessive junctional epithelium when treating periodontal defects, it is important to adapt Striate+™ closely to the treated tooth.
7. To prevent membrane displacement, Striate+™ may be fixed in place with sutures or pins, if required.
8. Use the previously created mucoperiosteal flap to close the wound over the membrane.
9. Complete wound closure is recommended but not essential. Excess tension to achieve wound closure may increase the risk of dehiscence and should be avoided.

### Post-operative care

1. Patients should be monitored closely in the initial post-operative period. The use of prophylactic antibiotics and oral antiseptics following surgery is recommended. Good oral hygiene is essential in the period following implantation and dental practitioners should provide additional guidance to patients on maintenance of oral hygiene post-treatment.
2. Post-operative symptoms may include swelling, pain or mild inflammation and dental practitioners should provide guidance to patients in appropriate symptom management.
3. Exposure of the Striate+™ membrane through wound dehiscence may occur and generally resolves spontaneously. Membrane removal is usually not required. In the event of membrane exposure, prophylactic treatment with antiseptic rinses to minimize the risk of bacterial contamination is recommended.

4. Allow sufficient time for bone regeneration before surgical re-entry.
5. Destructive parafunctional habits (bruxism, clenching), attrition or existing orthotic appliances may impede repair if in contact with the site of repair. In that case, dental practitioners should provide guidance to the patient in how to minimize damage to the treatment site.

### Precautions

Striate+™ should only be used by qualified dentists and oral surgeons trained in guided bone and tissue regeneration procedures.

Caution and close patient monitoring during the post-operative period may be required when using Striate+™ in patients with impacts on healing capacity due to:

- Uncontrolled metabolic disease (e.g. diabetes, thyroid disorders)
- Anti-coagulant / blood-thinning therapy
- Treatment with high doses of anti-inflammatory medications or bisphosphonates
- Connective tissue diseases
- Autoimmune diseases
- Radiotherapy
- Heavy smoking

### Adverse reactions

Complications that may be associated with the surgical procedure such as infection, dehiscence, membrane exposure, swelling, bleeding or pain should be discussed with the patient prior to the procedure.

Adverse reactions to porcine-derived collagen membranes in dental procedures are extremely rare, but immune reactions are possible.

Patients should be monitored closely in the initial post-operative period to identify and address any adverse reactions that may occur.

### Storage and handling

Store Striate+™ in its original packaging at controlled room temperature (15-25°C/59-77°F) in a dry place.



















### Presentation

Striate+™ is supplied in a double PETG/Tyvek tray contained in a labeled cardboard box.

Striate+™	
Ref	Size
OCG-152	15 x 20 mm
OCG-203	20 x 30 mm
OCG-304	30 x 40 mm

A paper IFU is available from your distributor, or you can make a request via the Orthocell website.

### Symbols used in labelling

	Use-by date
	Batch code
	Catalogue number
	Date of manufacture
	Importer
	Distributor
	Sterilized using irradiation
	Double sterile barrier system
	Do not sterilize
	Do not re-use
	Do not use if the blister packaging is damaged or opened.
	Keep dry
	Keep away from sunlight
	Temperature limit 15-25°C/59-77°F
	Contains biological material of animal origin
	Consult instructions for use
	Medical Device
	Unique device identifier



### Manufacturer

Orthocell Ltd, Building 191 Murdoch University,  
South Street, Murdoch WA 6150 Australia  
**orthocell.com | +61 8 9360 2888**

### Registration Holder:

S.I.N. IMPLANT SYSTEM LTDA  
Rua Soldado Ocimar Guimarães, 421- Vila Rio Branco.  
CEP: 03348-060 - São Paulo - SP - Brasil  
CNPJ : 04.298.106/0001-74  
Resp. TÉCNICO: Alessio Di Risio CREA-SP:  
5061207169  
Registro ANVISA 80108910132  
SERVIÇOS AO PROFISSIONAL  
SAC: 0800 770 8290

### SERVIÇOS AO PROFISSIONAL

SAC: 0800 770 8290

**www.sinimplantsystem.com.br**

E-mail: sin@sinimplantsystem.com