ortho·cell

Striate  TM

Guided Bone and Tissue Regeneration- Dental

advancing tissue
repair and regeneration

Striate+™ Product Information

Orthocell Striate+™

Bioabsorbable Implantable Collagen Membrane

Striate+™ is a type I collagen membrane which is sterile and suitable for use in guided tissue repair and bone regeneration in dental applications.

Striate+™ is ideal for use in surgical implantations as it is highly biocompatible with little antigenicity. The membrane enables the ingrowth and proliferation of cells and actively functions in normal physiological repair mechanisms. The membrane structure retains its strength as a tissue support when wet. The collagen membrane is a hydrophilic tissue support which conforms to the contours of the tissue surface. As the product is naturally derived, slight variations in the surface texture of the membrane are to be expected and these do not alter the functionality of the product. Each membrane is subject to careful inspection to assure high quality.

COMPOSITION

Striate+™ is a type I collagen medical device manufactured by a quality

controlled process. The raw material used for manufacture is approved by veterinarians to assure quality. The device has been manufactured and purified to remove reactive materials and to minimise the potential for immunologic reaction. The product is completely natural without the addition of additives or cross-linking to facilitate normal physiological processes of integration, resorption and remodelling.

Striate+™ has a bilayer structure with a porous and smooth layer. The porous (or rough) layer is to be applied contiguous to the tissue surface and promotes the ingrowth of cells. The smooth side facilitates smooth articulation within joints and cavities, preventing ingrowth of fibrous tissues.

Striate+™ is sterilised by γ -irradiation in sealed double-layer trays.

As a biologically derived product being predominantly comprised of collagen, Striate+™ is not anticipated to result in reciprocal interference during specific investigations or treatment.

INDICATIONS FOR USE

Striate+™ is intended as a tissue support and barrier membrane in guided bone and tissue regeneration procedures.

Striate+™ is recommended for:

- Simultaneous use of GBR-membrane (Striate+™) and implants.
- Augmentation around implants placed in immediate extraction sockets.
- Augmentation around implants placed in delayed extraction sockets.
- Localised ridge augmentation for later implantation; Alveolar ridge reconstruction for prosthetic treatment.
- Filling of bone defects after root resection, cystectomy, removal of retained teeth.
- Guided bone regeneration in dehiscence defects; and
- Guided tissue regeneration procedures in periodontal defects.

CONTRAINDICATIONS

Striate+™ must not be used if there is evidence of active infection at the treatment site. Striate+™ is

contraindicated for use in:

- Pregnant or lactating women.
- Immunocompromised or immunosuppressed patients.
- Patients with existing, untreated infections, including sinus or periapical pathology.
- Patients with prior malignancies of the head and neck treated with radiotherapy or antineoplastic / chemotherapy drugs.
- Patients being treated with medications affecting tissue metabolism (calcium channel blockers, epilepsy medications, cyclosporin).
- Patients with diseases affecting tissue metabolism (gingival fibromatosis or other fibrous dysplasias).
- Patients with diseases which may contribute to a poor healing response (gingival recession and other mucocutaneous disease, hypertension, uncontrolled diabetes, renal disease, collagen disorders [e.g. Ehlers-Danlos Syndrome]).
- Patients receiving treatment with high doses of NSAID medications or other medications which affect connective tissue metabolism (bisphosphonates). Maintained

usage at therapeutic plasma levels should be avoided for 3 months post-implant.

- Patients who have conditions causing bone or tissue loss (osteomalacia, osteoporosis, renal osteodystrophy).
- Patients with blood clotting disorders or anticoagulant / blood-thinning medication unless managed appropriately.
- Patients with clinically significant response to porcine materials or collagen.

DIRECTIONS FOR USE

General principles of surgical practice and aseptic technique must be adhered to when using Striate+™. All reasonable precautions to minimise the potential for contamination, including the control of saliva, must be taken.

The directions for use of Striate+™ are as follows:

1. Following the completion of surgical procedure(s), bone defects are filled as required with bone graft or other void-filling material (allograft, HA, natural bone substitute, e.g. Bio-Oss®).
2. Striate+™ is trimmed to the required size with sterile technique.
 - a. The membrane should significantly

overlap the walls of the defect to assure adequate enclosure and prevent soft-tissue invasion.

- b. The rough side of the membrane is placed facing the bone defect.
3. Striate+™ is applied over the defect and pressure applied until bleeding is adequately controlled and the membrane is uniformly wet with fluids (thus conforming to the defect void).
 4. Complete wound closure is recommended but not essential as excess tension may increase the risk of dehiscence.

POST-OPERATIVE CARE

Patients should be monitored in the initial post-operative period (4-6 weeks). Orthocell recommends the use of prophylactic antibiotics, oral antiseptics and/or education in post-operative oral hygiene in the weeks following surgery. Post-operative symptoms may include swelling, pain or mild inflammation. Surgeons are encouraged to instruct patients in appropriate symptom management and to make contact if concerned.

1. Surgical re-entry is not recommended for a period of at least 4 months to allow sufficient time for bone regeneration.

2. The membrane can be removed in the case of severe infection or adverse reaction. Inflamed or infected tissues should be excised with the residual membrane.
3. Prophylactic treatment with antiseptic rinses is highly recommended if the membrane becomes exposed due to wound dehiscence.

PRECAUTIONS

Striate+™ is intended for use only in dental and oral surgery by qualified dentists/oral surgeons.

Safety and efficacy in applications other than those indicated cannot be assured.

The long-term safety and effectiveness of the membrane in conjunction with bone-filling materials has not yet been established.

Use of the device in the presence of infection or inadequate oral hygiene may predispose the patient to an adverse reaction. Bone fillers or implants must be adequately localised or fixed in place (as appropriate) prior to use of Striate+™.

Safety and efficacy of repeat doses cannot be assured.

Safety and efficacy in patients exceeding the recommended age range of 18-65 years cannot be assured.

Destructive parafunctional habits (bruxism, clenching), attrition or existing orthotic appliances may impede repair if in contact with the site of repair. Dental practitioners should provide additional guidance to prevent damage to the repair.

Impeccable dental hygiene is a necessity in the period following implantation, smokers may need to modify their intake. It is recommended that practitioners ensure the likelihood of patient compliance to oral hygiene requirements.

Striate+™ has a smooth and a rough surface and the rough side is to be implanted in contact with the bony defect.

ADVERSE REACTIONS

Striate+™ is a porcine-derived collagen device. Although improbable due to decellularisation and removal of genetic material, persons highly sensitive to porcine material may react to the product. Humoral immunity against collagen is rare but can be detected

with a serologic test if concerns exist. General complications such as swelling, bleeding, dehiscence, pain etc. may be associated with the surgical procedure itself rather than the device and should be discussed with the patient prior to the procedure.

STORAGE AND HANDLING

Do not use the product if damaged or opened. The blister pack is tamper evident and cannot be opened without being visible.

The device is a single-use product. Re-use may result in severe adverse reactions, illness or infection. Any unused product must be discarded.

Use the product within expiry period stated on labels.

Striate+™ poses no known risks in all MR environments.

The product must be stored in its box (away from sources of light) at room temperature (15-25°C) in a dry place.

The device must be handled using sterile surgical technique.

APPEARANCE

Striate+™ is presented in a sterile double blister pack. Each package contains one (1) collagen membrane. The product codes for re-ordering are as follows.

Product No.	Size (mm)
OCG-152	15x20
OCG-203	20x30
OCG-304	30x40
OCG-405	40x50
OCG-050	50 (Ø/diameter)
OCG-070	70 (Ø/diameter)

SYMBOLS

Use-by date

Batch code

Catalogue number

Importer

Distributor

Sterilized using irradiation

Do not re-sterilize

Do not use if package is damaged

Double sterile barrier system

Keep away from sunlight

Keep dry

Temperature limit

Do not re-use

Consult instructions for use

SUPPLY INFORMATION

Authorized Representation in the European Union:

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Manufacturer:

Orthocell Ltd
Building 191 Murdoch
University
South Street Murdoch WA 6150
Australia
+61 8 9360 2888

CONCERNS / QUESTIONS

Contact Orthocell
info@orthocell.com.au

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MDSS GmbH